

CODE	DESCRIPTION	SIZE
T3603	Head and Neck	40x30cm
T3111	Head and Neck	50x40cm


MATERIALS: CUSHION IN MICROFIBER POLYESTER / PU (FULLY SEALED) AND INTERNAL POUCH (FULLY SEALED) WITH INSIDE THE TWO REAGENTS READY TO BE MIXED

BIOCOMPATIBILITY: ISO 10993-5:2009, 10993-12:2012, ISO 10993-10:2010

Do you want something custom-made?

We offer fully customized cushion solutions, tailored to meet each customer's specific needs. By working closely with our clients, we define the optimal cushion dimensions and personalize the fabric with their logo and brand colors.

Color: examples



*Colors shown are for reference only



ESSEBI MEDICAL SRL dba TOTIM[®]
Strada Campo del Fiume 84, 47896 Faetano - RSM
Email: info@totim.it

DESIGNED IN ITALY
Made in RSM

Copyright© Essebi Medical Srl Totim® Rev.02

OPHTHALMIC POSITIONING CUSHION

- IT MOULDS TO YOU -

**CUSTOMIZED CUSHION
THAT ADAPTS TO THE HEAD AND NECK
WITH A PERSONALIZED SHAPE**

Just hold to mold - easy to use - custom shaping

TRULY UNIQUE IN THE GLOBAL MARKET



MEMORY FOAM CUSHION FOR PATIENT POSITIONING DURING OPHTHALMIC SURGERY

TOTIM[®] Patient Cushion Positioning System is used for positioning adult and pediatric patients undergoing eye surgery in order to ensure more comfort to the patient and the surgeon

Cushion with external cover in microfiber polyester comfortable and absorbent water repellent non-woven material, markable, washable/easily sanitized - temperature controlled exothermic reaction - the short window of time required for solidification allows for the correct modeling - Single Patient Use

OCULAR SURFACE SURGERY

LIMBAL STEM CELL TRANSPLANTATION

CORNEAL SURGERY, INCLUDING FEMTOSECOND LASER SURGERY

COMPLEX CATARACT SURGERY OR IN HIGH-RISK PATIENTS

GLAUCOMA SURGERY

ANTERIOR SEGMENT SURGERY

POSTERIOR SEGMENT SURGERY

STRABISMUS SURGERY AND USE OF BOTULINUM TOXIN

OPHTHALMIC SURGERY IN CHILDREN UNDER GENERAL ANAESTHESIA

PATIENT HAS AN ALTERATION OF THE SPINE OR NECK: TRENDELENBURG POSITION

PATIENT CANNOT HAVE CHEST FLAT AND CANNOT EXTEND NECK: FACE-TO-FACE UPRIGHT SEATED POSITIONING

