


Are You Losing Precious Samples Because of Your Vials?




Three Features of the Shimadzu LabTotal Vial

- 1 Reduces the Adsorption of Basic Compounds on the Surface of the Glass Vial**
► This improves quantitative accuracy in the analysis of trace basic compounds.
- 2 Included Quality Certificate Signifies the High Quality of the Vial**
► This product can be used as a vial for high-sensitivity analysis in conventional LC and in LC/MS and GC/MS analysis.
- 3 Wide Mouth with Preset Cap and Septum Improves Ease of Use**
► This shortens the time needed for sample preparation and reduces human error.




Reduced Adsorption of Basic Substances

Adsorption of basic substances, a drawback of glass vials, has been reduced as much as possible. As a result, stable recovery rates can be obtained with no changes over time.



Mass Spec Quality Certificate Provided

This confirms that there was an absence of elution components from the vial in random inspections using LC/MS/MS and GC/MS/MS. Therefore, this product can be used with confidence, with no concern for ghost peaks originating from the vial.



Wide-Mouthed Vial, with Preset Septum and Cap

This product is designed so that the corners of the bottom of the vial interior can be reached by inserting a 100 μ L scale micropipette tip obliquely, which increases the ease of use when inserting and retrieving samples.

Also, marking labels, handy for note taking, are included with the vials.

P/N	Product Name
227-34001-01	LabTotal Vial Certified Kit (100 pcs)



• All prices indicated are excluding consumption tax.
• Installation and adjustment expenses will be charged separately.



Shimadzu Corporation
www.shimadzu.com/an/

Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

For Research Use Only. Not for use in diagnostic procedures.
The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.

© Shimadzu Corporation, 2015